

Defibrillators: Help to Save a Life

Early intervention of CPR and defibrillation from bystanders of cardiac arrest can often be the difference between life and death*. Ensure you are prepared in the event of cardiac arrest by exploring Banner's range of defibrillators and accessories.



Smarty Saver Semi-Automatic Defibrillator



Code	Description
1197933	Smarty Saver Semi-Automatic Defibrillator with Sturdy Defibrillator Case

Mediana A16 HeartOn Semi-Automatic AED



Code	Description
1170870	Mediana A16 HeartOn Semi-Automatic AED

Mediana A15 Battery Pack for HeartOn AED



Code	Description
1170869	Mediana A15 Battery Pack for HeartOn AED

Smarty Saver Fully Automatic Defibrillator



Code	Description
1197934	Smarty Saver Fully Automatic Defibrillator with Sturdy Defibrillator Case

Mediana A16 HeartOn Fully-Automatic AED



Code	Description
1170871	Mediana A16 HeartOn Fully-Automatic AED

Reliance Medical Mediana A15 HeartOn AED



Code	Description
1156061	Reliance Medical Mediana A15 HeartOn AED



CALL

1. Call 999



PUSH

2. Start CPR



SHOCK

3. Administer Defibrillator

How to Use a Defibrillator

- Step 1:** Press the green button to switch on the defibrillator and follow the instructions.
- Step 2:** Remove the person's clothing above the waist.
- Step 3:** Peel off the sticky pads and attach them to the person's bare skin. Put one pad on each side of the chest.
- Step 4:** Once you have attached the pads, stop CPR and do not touch the person. The defibrillator will then check the person's heart rhythm.
- Step 5:** The defibrillator will decide whether a shock is needed. If so, it will tell you to press the

'shock' button. An automatic defibrillator will shock the person without you needing to do anything. Do not touch the person while they're being shocked.

- Step 6:** The defibrillator will tell you when the shock has been given and whether you need to continue CPR.
- Step 7:** If the defibrillator tells you to continue to do CPR, continue with chest compressions until the person shows signs of life, or the defibrillator tells you to stop so it can analyse the heartbeat again.

* Judy O'Sullivan, British Heart Foundation Director of Innovation in Health Programs (2023)

* British Heart Foundation (2025)

Defibrillators Accessories



Products:

AED Alarmed Outdoor Wall Mountable Heated Metal Cabinet

- Thick carbon steel offers outdoor weather resistance
- Shatterproof transparent window for quick status checking
- User regulated temperature settings



Code	Description
1170876	AED Alarmed Outdoor Wall Mountable Heated Metal Cabinet Yellow

Mediana AED Metal Wall Cabinet

- Toughened durable glass door
- Alarmed cabinet sounds when opened
- Lockable and includes keys



Code	Description
7089887	Mediana AED Metal Wall Cabinet with Glass Door and Alarm Lockable Large

Blue Dot AED Emergency Response Kit

- Ideal kit to pair with defibrillators
- Contains alcohol free cleansing wipes, resuscitation face shield and a pair of gloves
- Also contains sterile gauze dressing pads, disposable razor and a pair of scissors



Code	Description
1465039	Blue Dot AED Emergency Response Kit

For further information please contact your Banner account manager or email contact.banner@BannerUK.com

